




BindPro™

Protein Removal for Drug Binding/Screening and Metabolomics

- Linearly scaleable, unlike ultrafiltration
- Suitable for use with surfactants, chaotropes, water-soluble analytes
- Fast process, less than 15 minutes from application to separation
- Applicable for drug binding/screening, metabolomics and protein recovery

BindPro™ is a polymeric protein removal suspension reagent. It is designed as an alternative to ultrafiltration for applications that require a more versatile or scaleable format. BindPro™ also can be used in lieu of solvents for drug binding studies, especially useful for analytes that are water soluble. Consequently, BindPro™ has applications in a range of drug binding, screening and metabolomic investigations. If desired, proteins can be recovered from BindPro™ under moderately alkaline conditions.



| Protein | BindPro™: Sample | Removal |
|--------------------------|------------------|---------|
| BSA, PBS @ 30 mg/ml | 1 : 1 | >99% |
| BSA, 1%SDS @ 30 mg/ml | 1 : 1 | >99% |
| BSA, 3M GuSCN @ 30 mg/ml | 1 : 1 | >99% |
| Human Serum | 2 : 1 | >99% |



BIOTECH SUPPORT GROUP

| Product | Size | # of Samples & Sample Size* | Item No. | Price |
|-----------------|-------|-----------------------------|----------|-------|
| BindPro™ | 15 ml | 75, 100µl Serum Samples | BP355-15 | \$325 |
| BindPro™ | 50 ml | 250, 100µl Serum Samples | BP355-50 | \$680 |

| Items Required | 10 Prep | 50 Prep | Reagent |
|-------------------------------|----------|-----------|-----------------|
| BindPro™ | 0.5 gram | 2.5 grams | Supplied |
| Binding Buffer HVBB, PH 6.0 | 15 ml | 125 ml | Supplied |
| Wash Buffer HVWB, PH 7.0 | 15 ml | 125 ml | Supplied |
| Elution Buffer HVEB, PH 9.8 | 15 ml | 125 ml | Supplied |
| SpinX Centrifuge tube filters | 10 | 50 | Supplied |

PROTOCOL

1. Resuspend BindPro™ by shaking well prior to use.
2. Add 2ml of BindPro™ to 1 ml of the sample (2:1 volume ratio). Use wide bore pipette tips.
3. Gently mix by inversion for 10 minutes at room temperature.
4. Centrifuge sample at 10,000 x g for 5 minutes or microfuge at 16,000 x g for 5 minutes.
5. Retain the precipitate, which contains bound proteins and is ready for further processing.

Note: To elute the bound proteins, a high pH buffer can be employed (pH >9.0)

References

Wolbachia Lipoprotein Stimulates Innate and Adaptive Immunity through Toll-like Receptors 2 and 6 to Induce Disease Manifestations of Filariasis Joseph D. Turner J. Biol. Chem. 2009 284: 22364-22378.

CONTACT US

We welcome your questions and comments regarding our products.

Call 732-274-2866, 800-935-0628 (North America) Mon – Fri 9am-6pm EST.

Fax 732-274-2899

Email sales@biotechsupportgroup.com

Mail 1 Deer Park Drive, Suite P, Monmouth JCT, NJ 08852